

THERMOCOUPLE SENSOR



Thermocouple Sensor

Part number 10B-Class 1-T-5-CAP. Class 1 type T exposed junction sensor, 1/0.3mm ptfе flat pair leads 5m long. Sensor end cap reduces water ingress.



Thermocouple Entry Gland

Part number SVMS 264. Small entry gland. Stainless steel and brass construction, 12 thermocouples with pressure port; 1/4" BSP fitting.
Part number SVMS 266. Large entry gland. Stainless steel and brass construction, 30 thermocouples with pressure port; 3/4" BSP fitting.

One of a range of products from Anville Instruments for sterilizer and washer-disinfector validation to meet standards HTM0101, HTM2010, HTM2030, EN285 & EN584.



Anville Instruments Ltd
Bramble Orchard, 46a Folly Lane North,
Farnham, Surrey, GU9 0HX.

e-mail: sales@anvilleinstruments.com

Main No. +44 (0) 1252 351 030.

Fax: +44 (0) 1252 323 492.

